

## 300 kN DUAL COLUMN FLOOR MODEL

CP128610

CP128610 systems perform tensile, compression, bend, creep, and cyclic tests on all raw materials and finished goods. These testing instruments are engineered for precision, built for durability, and offer flexibility for changing requirements. They are designed with features that increase testing efficiency and improve the testing experience for the operator.

CP128610 floor model systems are robust, heavy-duty frames commonly used for testing high-strength metals and alloys, advanced composites, aerospace and automotive structures, bolts, fasteners, and plate steels. A key advantage of the 5980 floor models is the stiffness of the frame that is especially important when testing high-strength materials such as aerospace composites, metal alloys, and crystalline polymers. Pre-loaded bearings and precision ball screws, a thick crosshead and base beam, and low-stretch drive belts contribute to better performance by producing more accurate modulus and strain values, and minimizing the energy stored during a test.

Easier access to the test area is a very important benefit for the operator. When testing large metal specimens, fasteners, or composite parts, operators often struggle to reach the test area of large frames in order to load a specimen - especially when the grips or fixtures are large. These systems are designed with a lower base height that allows the operator to stand closer to the grips and fixtures, reducing the need to hold heavy test specimens away from the body.

### FEATURES

- Meets or exceeds requirements of all national and international standards; namely ISO, ASTM, BS, DIN, EN, and AFNOR
- Thousands of accessories to meet test requirements in almost any application or industry: biomedical, automotive, electronics, plastics, metals, composites, elastomers, aerospace, textiles, and many more
- Supported by the largest global Service organization in the industry; delivering high-quality calibrations, training, preventative maintenance, and technical support
- Productivity Panel with Live Display, Soft keys and Specimen Protect for enhanced usability and productivity
- Power and electrical components for improved shock load performance
- Shock mounted electrical components

### BLUEHILL® UNIVERSAL AND INSTRON® CONNECT

Designed from the ground up for touch, Instron's static testing software, Bluehill Universal, is easy-to-use, increases testing efficiency, and contains modular features that enable users to run the most complex of tests.

With ISO 9001 accreditation, our goal is to provide the best ownership experience by delivering the highest quality products, expert support, and world-class service. Instron Connect provides users with a powerful communication platform via a secure connection between the Instron system at your facility and Instron's global technical support engineers. With Instron Connect, users receive faster remote technical support, reduce risk with schedule verification and preventive maintenance reminders, and are effortlessly able to keep up to date with the latest software features.



## SPECIFICATIONS

Force Capacity <sup>1</sup>	kN	300
	lbf	67400
Vertical Test Space <sup>3</sup>	mm	1830
	in	72
Horizontal Test Space <sup>4</sup>	mm	575
	in	22.6
Testing Speed Range Min-Max (Return)	mm/min	0.00005 - 400 (480)
	in/mm	0.000002 - 15.7 (18.9)
Position Control Resolution	nm	9
	µin	0.33
Frame Axial Stiffness	kN/mm	350
	lbin	1998600
Maximum Force at Full Speed	kN	150
	lbf	33720
Maximum Speed at Full Force	mm/min	200
	in/min	8
Height	cm	278
	in	110
Width <sup>5</sup>	cm	113
	in	44
Depth	cm	78
	in	31
Weight	kg	955
	lbs	2100
Maximum Power Requirement <sup>7</sup>	VA	350

## Notes:

1. Meets or exceeds ASTM E4, BS 1610, DIN 51221, ISO 7500/1, EN 10002-2, JIS B7721, JIS B7733, and AFNOR A03-501 standards. Instron recommends that systems are verified on-site at the time of installation as required by ASTM E4 (par. 20.3) and ISO 7500-1 section 9) standards.
2. All systems conform to all relevant European standards and carry a CE mark.
3. Vertical test space on all systems is the distance from the top surface of the base platen to the bottom surface of the moving crosshead, excluding load cell grips and fixtures.
4. Horizontal test space on dual column systems is the distance between the inside edge of the columns.
5. This is the system footprint width. The Operator Dashboard monitor may add 300 mm (12 in) to the overall width of the frame.
6. These specifications were developed in accordance with Instron's standard procedures and are subject to change without notice.
7. CP128610 systems are single phase

## Additional Specifications

**Force Measurement Accuracy:**

± 0.4% of reading down to 1/100 of load cell capacity with 2525, 2530 or 2580 Series load cells; ± 0.5% of reading down to 1/1000 of load cell capacity with 2580 Series load cells (with Advanced Performance Option); ± 0.5% of reading down to 1/500 of load cell capacity with 2580 Series load cells; ± 0.5% of reading to 1/250 of load cell capacity with 2525 or 2530 Series load cells. <sup>1</sup>

**Displacement Measurement Accuracy:**

±0.01 mm or 0.05% of displacement (whichever is greater)

**Strain Measurement Accuracy:**

Meets or surpasses the following standards: ASTM E83, ISO 9513, and EN 10002-4.

**Testing Speed Accuracy:**

(Zero or constant load): ±0.1% of set speed

**Data Acquisition Rate at the PC:**

Up to 2.5 kHz (Advanced option) simultaneous on force, displacement, and strain channels, 1 kHz (Standard)

## Facility Requirements and Operating Environment

**Single Phase Voltage:**

200/208, 220, 230, 240 VAC ±10%, 47 to 63 Hz. Power supply must be free of spikes, surges or sags exceeding 10% of the average voltage.

**Three Phase Voltage:**

208, 240, 400 or 480 VAC ±10%; 47 to 63 Hz. Power supply must be free of spikes, surges or sags exceeding 10% of the average voltage.

**Operating Temperature:**

+10 to +38 °C (+50 to +100 °F)

**Storage Temperature:**

-40 to +66 °C (-40 to +150 °F)

**Humidity Range:**

+10 to +90%, non-condensing

**Atmosphere:**

Designed for use under normal laboratory conditions. Protective measures may be required if excessive dust, corrosive fumes, electromagnetic field, or hazardous conditions are encountered.

[www.instron.com](http://www.instron.com)



Worldwide Headquarters  
825 University Ave, Norwood, MA 02062-2643, USA  
Tel: +1 800 564 8378 or +1 781 575 5000

European Headquarters  
Coronation Road, High Wycombe, Bucks HP12 3SY, UK  
Tel: +44 1494 464646